### **LINISYN** A division of GE Healthcare

# Your Partner for multi-brand Probe Repair

Unisyn<sup>\*</sup>, a division of GE Healthcare, is your trusted provider for ultrasound probe solutions.

Unisyn supports your ultrasound equipment by providing comprehensive, convenient, and economic probe repair and exchange<sup>1</sup>.



- In-depth FirstCall\* probe testing services
- Extensive probe repair capabilities
- A variety of options for high quality probe repair and exchange<sup>1</sup>
- Convenient access to probe loaners
- Rapid probe evaluation and repair turnaround



## FirstCall<sup>\*</sup> probe testing, integrated into our process

Proprietary to Unisyn<sup>\*</sup>, our FirstCall probe testing device provides a detailed analysis of probe performance. FirstCall testing is integrated into our inspection, diagnosis, and repair processes. We use FirstCall technology to verify each customer repair, as well as probes we add to our inventory, each step of the way.

#### FirstCall probe diagnosis

During the inspection process, FirstCall can properly identify common ultrasound probe failures, such as:

- The total number and exact location of dead or weak elements within the array
- Changes in the acoustic performance of the array
- Probe lens damage and delamination
- Broken wires within the probe cable, flex circuit, and probe head
- Defective electronics within the probe connector.

#### Inventory you can rely on

We stock a large quantity of loaner probes across several manufacturers and we check the quality of our probe inventory using our FirstCall technology. We know it is critical to deliver probes to you as quickly as possible so you can get your ultrasound equipment up and running.

#### Thorough probe inspection and evaluations

All of the probes that come through our facility go through an in-depth inspection process to diagnose and evaluate clinical performance as well as assess the cosmetic condition of the probe. Our extensive four-point inspection process provides you with a detailed analysis of the health of your probe.

You will receive a detailed four-point evaluation summarizing the condition of your probe including:

- Visual inspection of the probe's cosmetic condition
- Electrical leakage testing for the probe's structural and electrical integrity
- FirstCall testing that measures the probe's electrical and acoustic properties
- System and console scanning to check the overall imaging performance (available for most probe models).



#### Contact us: 0845 850 3392

Unisyn – GE Healthcare 71 Great North Road Hatfield - Herts AL9 5EN



©2013 General Electric Company – All rights reserved.

#### Data subject to change.

Marketing Communications GE Medical Systems Société en Commandite Simple au capital de 85.418.040 Euros 283 rue de la Minière – 78533 Buc Cedex France RCS Versailles B 315 013 359

A General Electric company, doing business as GE Healthcare

GE and GE Monogram are trademarks of General Electric Company.

\* Trademarks of General Electric Company. 1 Probe exchange available for non-repairable GE probes only.

